




Quantitative Reference Cultures USP 61/62, Growth Promotion and Inhibitory Testing



Innovative Controls to Support Data Quality

Traceable

Manufactured from USP specified pure organism strains using media and procedures conforming to USP specifications.

Easy To Use

Supplied with sterile hydration fluid, simply add a single sample to the 1.2 mL buffer and swirl. Resulting inoculant activity is less than 100 CFU/0.1 mL.

Certified

The activity of each lot of organisms is analytically certified by spiral count on USP specified media. We analyze at least 20 samples from each production lot for activity confidence.

Economical

Each organism strain is supplied in a kit of 10 individual samples and 10 vials of sterile hydration fluid at a value price point.

Quantitative Reference Cultures

USP 61/62, Growth Promotion and Inhibitory Testing

Organism	Strain Identification	Activity	10 Vials Catalog Number/Price
S. aureus	ATCC6538/NCTC10788/WDCM0032	~1000 CFU	61-0032/\$225.00
P. aeruginosa	ATCC9027/NCTC12924/WDCM00026	~1000 CFU	61-00026/\$225.00
B. subtilis	ATCC6633/NCTC10400/WDCM0003	~1000 CFU	61-0003/\$225.00
C. albicans	ATCC10231/NCPF3179/WDCM00054	~1000 CFU	61-00054/\$225.00
A. brasiliensis	ATCC16404/NCPF2275/WDCM00053	~1000 CFU	61-00053/\$225.00
E. coli	ATCC8739/NCTC12933/WDCM0012	~1000 CFU	61-0012/\$225.00
S. enterica	ATCC14028/NCTC12023/WDCM00031	~1000 CFU	61-00031/\$225.00
S. enterica abony	NCTC6017/WDCM00029	~1000 CFU	61-00029/\$225.00